


Download APP Planmeca ClinicPatrol™



Download APP

- **Scan QR-koden**

Open camera to scan QR code and use link



(QR-Code is also found in the user manual)


- **Download via link**

The APP can also be downloaded directly via this link:

www.bwt-group.net/bwt/download/bwt-RO-Systems.html


Install APP

Click 'Install'




Android devices

If using an Android device the APP is now ready




iOS-devices


- Go to settings on your device and choose 'General'




- Choose 'VPN & Device Management'




- Choose 'BWT Aktiengesellschaft'




- Click 'Trust BWT Aktiengesellschaft'



- Click 'Trust' again



The app is now installed on your device and ready to use



Recommended settings for pressure tanks

Tank size	2L.	8L	25L
Pre pressure for tank	1 Bar	1 Bar	1 Bar
Hysteresis	1,6 Bar	1,6 Bar	1,6 Bar
Stop pressure	3,0 Bar	3,0 Bar	3,0 Bar